

Quantity original	Bezeichnung Hydraulik - und Elektrikbauteile	Description Hydraulic & Electrical components	Siemens No.	Bezeichn.	Riedel No.	Supplier CSML	Supplier RIEDEL
1	Manometer	Pressure Gauge	10161725		4541038529	x	
1	Ausdehnungsgefass	Expansion vessel	8905304		4541038531	x	
1	Fuell-Entleerungshahn	Valve Drain	3056884		4541018514	x	
1	Kugelabsperrentil	Ball Valve	3054574		4525058504	x	
1	PTC-Heizung 30W	PTC crankcase heater 30 W	n/a		4541028520		x
1	Sekundär Pumpe CHI 4-50 incl. Kabel, Stecker und Verschraubung	Secondary Pump CHI 4-50 incl. Cable, plug and union (P30 Split)	8905379	M6	4541038537	x	
1	Primär Pumpe CHI 4-50 incl. Kabel, Stecker und Verschraubung	Primary Pump CHI 4-50 incl. Cable, plug and union	4807694	M3	4541028501	x	
1	Spirovent	Air purge	4820952		4541028508	x	
3	Temperatursensor	Temperature Sensor	7395432		4541038503	x	
1	Baugruppe Schaltraum WKS P30	Control Unit WKS P30	10189166		4541038538	x	
1	Regelplatine WKS	Control circuit board WCS	7395424	N1	4541038502	x	
1	Rueckkuehler-Ventilator	Cooler Fan	8905361		4541038539	x	
1	Rueckkuehler 2NK6 725-1A	Outdoor unit 2NK6 725-1A	8905346		2NK67251A	x	
1	Permanententluefter	Permanent Air purge valve	3056751		4541018504	x	
1	Kuehlwasserregler	Water regulating valve	8905387		4541038541	x	
1	Uberstroemventil	Pressure relief valve	8905395		4541038527	x	
3	Entstoerglied 3RT1916	Electr. noise suppressor	n/a		4541024109		x
1	Entstoerglied 3RT1926	Electr. noise suppressor	n/a		4541024115		x
1	Hilfsschalterblock	Contacteur auxiliary	n/a		4541024118		x
3	Hilfsschalter 3RV1901	Contacteur auxiliary	n/a		4541024122		x
3	Leistungsschalter	Circuit Breaker	3072410	F3,F10,F20	4541024120	x	x
3	Schuetz 3RT1015	Contacteur	3072378	K1,K3,K4,K5	4541024107	x	x
1	Schuetz 3RT1025	Contacteur	3072386	K2	4541024108	x	x
1	Sicherungsautomat D3	Circuit breaker (Fuse)	4694386	F8	4541024114	x	x
1	Sicherungsautomat	Circuit breaker (Fuse)	7743198	F9	4541024138	x	x
1	Transformator	Transformer	3073442	T1	4541034137	x	x
1	Schuetz 3RT1015	Contacteur	8682259	K20	4518034652	x	x
1	Entstoerbauteil fuer Schuetz 3RT 1916	Electr. noise suppressor for contactor	n/a		4541034123		x
1	Hilfsschalter	Contacteur auxiliary	n/a		4518034658		x
1	Kapillarrohrrregler	Capillary controller	4015988		4541038556	x	
	Kuehlkreislauf Komponenten	Refrigeration Circuit Components					
1	Heissgasbypassventil	Hotgas bypass valve	n/a		4518554106		x
1	Expansionsventil mit Isolierung	Expansion Valve with insulation	n/a		4541028514		x
1	Filterrockner	Filter drier	n/a		4541028515		x
1	Platten-WT	Plate heat exchanger	n/a		4541028516		x
1	Pressostat HD	Pressure switch HP	n/a		4541028517		x
1	Pressostat ND	Pressure switch LP	n/a		4541028519		x
1	Platten-WT Primär	Plate heat exchanger primary	n/a		4541034118		x
1	Rueckschlagventil NRV12S	Check valve Brass refrig.	n/a		4541028521		x
1	Rueckschlagventil NRV10S	Check valve Brass refrig.	n/a		4541028522		x
1	Schauglas SGN 12S	Sight Glass	n/a		4541028523		x
10	Schraderventileinsatz	Schrader valve insert	n/a		4541038535		x
1	Verdichter MTZ 50-4V	Compressor	n/a		4541028525		x
1	Verdichter MTZ 18-4V	Compressor	n/a		4541028526		x
1	Verfluessiger W/L	Condensor W/L	n/a		4541022801		x